



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz weighing glass with NS cap, Ø 48 x H. 70 mm, ground joint NS 50/12, high form**

**€196.68  
plus VAT &  
Shipping**

## Product Images



## Description

---

### Weighing glasses with NS cap

With replaceable NS cap, with inscription field.

Dimensions (mm) and contents (ml) permanently burned in.

## Additional Information

---

|                        |                     |
|------------------------|---------------------|
| No.                    | LZ-1538             |
| Manufacturer (Brand)   | Lenz                |
| old neoLab article no. | B-3541              |
| EAN                    | 4058072312770       |
| Transport temperature  | Room temperature    |
| diameter               | 48 mm               |
| Color                  | transparent (clear) |
| Material               | Borosilicate glass  |
| Volume                 | 110 ml              |
| Form                   | high shape          |
| Height                 | 70 mm               |
| Type closure/ lid      | NS cap              |

