



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz weighing glass with NS cap, Ø 75 x H. 95 mm, ground joint NS 80/12, high form

**€63.06
plus VAT &
Shipping**

Product Images



Description

Weighing glasses with NS cap

With replaceable NS cap, with inscription field.

Dimensions (mm) and contents (ml) permanently burned in.

Additional Information

No.	LZ-1539
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3542
EAN	4058072312787
Transport temperature	Room temperature
diameter	75 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	350 ml
Form	high shape
Height	95 mm
Type closure/ lid	NS cap

