



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Disposable weighing dishes white, 41 x 41 x 8 mm, 500 pcs./pack

€34.60
plus VAT &
Shipping

Product Images



Description

Flexible plastic weighing dishes made of blue plastic, for liquids and solids. With light-colored powder, it is easy to see whether there are still residues in the bowl.

Material: Polystyrene, temperature-resistant up to 82°C

Antistatic

Additional Information

| | |
|---------------------------------|--------------------------------|
| No. | 1-1124 |
| Manufacturer (Brand) | neolab |
| EAN | 4058072031954 |
| Color | White |
| Material | Polystyrene (PS) |
| Volume | 7 ml |
| Wide | 41 mm |
| Form | rectangular |
| Height | 8 mm |
| Length | 41 mm |
| Depth | 41 mm |
| TBST MAX | 82 °C |
| Bowls, pots and jars properties | stackable for 1xuse antistatic |
| Type edge | straight edge |
| Type bowl, dish, crucible | Weighing dishes |

