



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Disposable weighing dishes white, 89 x 89 x 25 mm, 500 pcs./pack

€53.20
plus VAT &
Shipping

Product Images



Description

Flexible plastic weighing dishes made of blue plastic, for liquids and solids. With light-colored powder, it is easy to see whether there are still residues in the bowl.

Material: Polystyrene, temperature-resistant up to 82°C

Antistatic

Additional Information

No.	1-1125
Manufacturer (Brand)	neoLab
EAN	4058072031961
Color	White
Material	Polystyrene (PS)
Volume	100 ml
Wide	89 mm
Form	rectangular
Height	25 mm
Length	89 mm
Depth	89 mm
TBST MAX	82 °C
Bowls, pots and jars properties	stackable for 1xuse antistatic
Type edge	straight edge
Type bowl, dish, crucible	Weighing dishes

