



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Weighing dishes diamond shape,
black, 25 ml, 71x46x14 mm, 500 pcs./pack**

**€66.30
plus VAT &
Shipping**

Product Images



Description

Flexible, antistatic polystyrene weighing boats in diamond shape facilitate pouring. The black color provides better contrast with light-colored weighing material.

PS material is resistant to weak acids, aqueous solutions, alcohol and alkalis. Temperature resistant up to 70°C.

Additional Information

No.	1-1375
Manufacturer (Brand)	neolab
EAN	4058072032807
Color	Black
Material	Polystyrene (PS)
Volume	25 ml
Wide	46 mm
Form	diamond-shaped
Height	14 mm
Length	71 mm
Depth	71 mm
TBST MAX	70 °C
Bowls, pots and jars properties	stackable/connectable for 1xuse
Type edge	straight edge
Type bowl, dish, crucible	Weighing dishes

