



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Weighing dishes diamond shape,
black, 100 ml, 119x90x19 mm, 500 pcs./pack**

**€112.10
plus VAT &
Shipping**

Product Images



Description

Flexible, antistatic polystyrene weighing boats in diamond shape facilitate pouring. The black color provides better contrast with light-colored weighing material.

PS material is resistant to weak acids, aqueous solutions, alcohol and alkalis. Temperature resistant up to 70°C.

Additional Information

No.	1-1377
Manufacturer (Brand)	neolab
EAN	4058072032821
Color	Black
Material	Polystyrene (PS)
Volume	100 ml
Wide	90 mm
Form	diamond-shaped
Height	19 mm
Length	119 mm
Depth	119 mm
TBST MAX	70 °C
Bowls, pots and jars properties	stackable/connectable for 1xuse
Type edge	straight edge
Type bowl, dish, crucible	Weighing dishes

