



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Disposable Weighing Dishes (PS), w. spout, 56 x 41 x 8 mm, 250 pcs./pack

€33.90
plus VAT &
Shipping

Product Images



Description

Drop-shaped disposable weighing trays made of PS.

Stable design, therefore also suitable for liquids.

Temperature resistant up to +93°C.

Antistatic.

One side is provided with a spout.

Additional Information

No.	1-1136
Manufacturer (Brand)	neoLab
EAN	4058072032012
Color	transparent (milky)
Material	Polystyrene (PS)
Volume	20 ml
Wide	44.5 mm
Form	drop-shaped
Height	7.93 mm
Length	57.2 mm
Depth	57.15 mm
TBST MAX	93 °C
Bowls, pots and jars properties	stackable/assemblable with spout for 1xuse antistatic
Type edge	straight edge
Type bowl, dish, crucible	Weighing dishes

