



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® autoclave tray, 60 x 30 x 26 cm, with 180° lid opening

**€649.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Container made of stainless steel (X5CrNi1810) for the transport of e.g. contaminated medical material.

Airtight lid with silicone seal, hermetically closed.

Stainless steel carrying handles with silicone coating.

Holder for identification labels.

Airtight and leak-proof closure.

Stackable.

Autoclavable (autoclave separately).

Riveted joints, no material changes due to welding.

Polished outer wall facilitates cleaning and maintenance.

Rounded corners for increased safety.

## Additional Information

---

No.	7-8186
Manufacturer (Brand)	neoLabLine
EAN	4058072168155
Transport temperature	Room temperature
Color	Silver Blue
Material	Stainless steel
Volume	37 ml
Wide	300 mm
Width/inside	260 mm
Form	rectangular
Height	260 mm
Height inside	260.0 mm
Length	600 mm
Length inside	560 mm
Bowls and tubs features	autoclavable with closure/lid
Type tub/ box	Autoclave tray

