




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® transport container stainless steel, leak-proof, 60x30x26 cm

**€1,066.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rigid transport container completely made of stainless steel (X5CrNi1810) for the transport of contaminated medical material.

Complies with EN ISO 11607-1 and EN 868-8, conforms to 93/42/EEC - Class I.

Airtight lid with silicone seal, hermetically closed.

Stainless steel carrying handles with silicone coating.

2 holders for aluminum or paper identification tags.

Tamper evident closure possible with holder for security seal.

Airtight and leak proof closure.

Stackable.

Autoclavable (autoclave separately).

Riveted joints, no material changes due to welding.

Polished outer wall facilitates cleaning and maintenance.

Rounded corners for increased safety.

Outer dimensions: 600 x 300 x 260 mm (L x W x H)

**Usable inner surface:** 560 x 260 x 250 mm (L x W x H)

## Additional Information

No.	7-8180
Manufacturer (Brand)	neoLab
EAN	4058072042233
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	46.8 l
Wide	300 mm
Width/inside	260 mm
Form	angular
Height	260 mm
Height inside	250.0 mm
Length	600 mm
Depth	600 mm
Bowls and tubs features	with special chemical resistance with closure/lid with handle suitable for foodstuffs
Length inside	560 mm
Type tub/ box	Transport tray/box

