



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Usbeck drip tray stainless steel, 310 x 125 x 30 mm

€101.30
plus VAT &
Shipping

Product Images



Description

The ideal stand for placing hot items on the laboratory bench.

Also for frozen containers that are to thaw slowly on the lab bench.

The liquid is collected in the tray and can be emptied if required.

The perforated plate can be removed.

Additional Information

No.	UB-0394
Manufacturer (Brand)	Usbeck
old neoLab article no.	1-1190
EAN	4058072401955
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	1.2 l
Wide	125 mm
Form	angular
Height	30 mm
Length	310 mm
Depth	310 mm
Bowls and tubs features	with special chemical resistance suitable for foodstuffs
Type tub/ box	Drip tray

