



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Usbeck drip tray stainless steel, 320 x 200 x 27 mm

**€111.65**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The ideal stand for placing hot items on the laboratory bench.

Also for frozen containers that are to thaw slowly on the lab bench.

The liquid is collected in the tray and can be emptied if required.

The perforated plate can be removed.

## Additional Information

---

No.	UB-0395
Manufacturer (Brand)	Usbeck
old neoLab article no.	1-1191
EAN	4058072401962
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	1.7 l
Wide	200 mm
Form	angular
Height	27 mm
Length	320 mm
Depth	320 mm
Bowls and tubs features	with special chemical resistance suitable for foodstuffs
Type tub/ box	Drip tray

