




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Ice bucket with lid, approx. 2.5 l, PU, black

**€144.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Ice bucket made of foamed PU, manufactured in one piece.

- thick-walled with good insulating properties
- complete with lid
- suitable for ice and dry ice as well as for liquid nitrogen
- Proven US design

Temperature resistant from -196°C to +100°C

## Additional Information

---

No.	2-6002
Manufacturer (Brand)	neoLab
EAN	4058072004842
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	268 mm
Color	Black
Material	Polyurethane (PU)
Volume	2.5 l
sterile	No
Height	179 mm
Refrigerant medium	Ice liquid nitrogen cold mixture dry ice
TBST MAX	100 °C
TBST MIN	-196 °C
Bowls and tubs features	with closure/lid insulating grid insert available as accessory

