



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Ice bucket with lid, approx. 2.5 l, PU,
red**

**€138.60
plus VAT &
Shipping**

Product Images



Description

Ice bucket made of foamed PU, manufactured in one piece.

- thick-walled with good insulating properties
- complete with lid
- suitable for ice and dry ice as well as for liquid nitrogen
- Proven US design

Temperature resistant from -196°C to +100°C

Additional Information

No.	2-6004
Manufacturer (Brand)	neoLab
EAN	4058072004859
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	268 mm
Color	red
Material	Polyurethane (PU)
Volume	2.5 l
sterile	No
Height	179 mm
Bowls and tubs features	with closure/lid insulating grid insert available as accessory
Refrigerant medium	Ice liquid nitrogen cold mixture dry ice
TBST MAX	100 °C
TBST MIN	-196 °C

