




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Lightweight ice tray without lid, 1 l, PVC


Special Price

€60.00 was

€75.00

**plus VAT &
Shipping**

Product Images



Description

The neoLab ice bins are characterized by their light weight and handiness.

They are made of PVC, which makes the containers much lighter than many other ice containers on the market.

With a temperature resistance of -196°C to 93°C, the neoLab ice containers are versatile.

The ice containers do not deform and are ideal for use with liquid nitrogen, dry ice, alcohol or salt solutions.

Additional Information

No.	2-4991
Manufacturer (Brand)	neoLabLine
EAN	4058072122201
Transport temperature	Room temperature
Color	yellow
Material	Polyvinyl chloride (PVC. non-specific)
Volume	1 l
Wide	145 mm
Form	rectangular
Height	80 mm
Length	183 mm
Depth	183 mm
Refrigerant medium	Ice Dry ice Liquid nitrogen
TBST MAX	93 °C
TBST MIN	-196 °C
Bowls and tubs features	insulating

