




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Lightweight ice bucket with lid, vol. 4 l, PVC

Special Price

€84.00 was

€105.00

**plus VAT &
Shipping**

Product Images



Description

The neoLab ice bins are characterized by their light weight and handiness.

They are made of PVC, which makes the containers much lighter than many other ice containers on the market. With a temperature resistance of -196°C to 93°C, the neoLab ice containers are versatile.

The red ice containers do not deform and are ideal for use with liquid nitrogen, dry ice, alcohol or salt solutions.

Additional Information

No.	2-4992
Manufacturer (Brand)	neoLabLine
EAN	4058072080686
Transport temperature	Room temperature
diameter	268 mm
Color	red
Material	Polyvinyl chloride (PVC. non-specific)
Volume	4 l
Form	rectangular
Height	227 mm
Refrigerant medium	Ice Dry ice Liquid nitrogen
TBST MAX	93 °C
TBST MIN	-196 °C
Bowls and tubs features	with closure/lid insulating

