



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lightweight ice tray without lid, Vol. 4 l, PVC

Special Price
€79.00 was
€98.10
plus VAT &
Shipping

Product Images



Description

The neoLab ice bins are characterized by their light weight and handiness.

They are made of PVC, which makes the containers much lighter than many other ice containers on the market.

With a temperature resistance of -196°C to 93°C, the neoLab ice containers are versatile.

The green ice containers do not deform and are ideal for use with liquid nitrogen, dry ice, alcohol or salt solutions.

Dimensions:

Length inside bottom: 235 mm

Length inside top: 260 mm

Width inside bottom: 150 mm

Width inside top: 170 mm

Additional Information

No.	2-4993
Manufacturer (Brand)	neoLabLine
EAN	4058072080693
Transport temperature	Room temperature
Color	green
Material	Polyvinyl chloride (PVC. non-specific)
Volume	4 l
Wide	250 mm
Form	rectangular
Height	127 mm
Length	350 mm
Depth	350 mm
Refrigerant medium	Ice Dry ice Liquid nitrogen
TBST MAX	93 °C
TBST MIN	-196 °C
Bowls and tubs features	insulating

