




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Lightweight ice tray without lid, vol. 9 l, PVC

Special Price

**€113.00** was

**€142.00**

**plus VAT &  
Shipping**

### Product Images



## Description

---

**The neoLab ice bins are characterized by their light weight and handiness.**

They are made of PVC, which makes the containers much lighter than many other ice containers on the market.

With a temperature resistance of -196°C to 93°C, the neoLab ice containers are versatile.

The yellow ice containers do not deform and are ideal for use with liquid nitrogen, dry ice, alcohol or salt solutions.

## Additional Information

---

No.	2-4994
Manufacturer (Brand)	neoLabLine
EAN	4058072080709
Transport temperature	Room temperature
Color	yellow
Material	Polyvinyl chloride (PVC. non-specific)
Volume	9 l
Wide	355 mm
Form	rectangular
Height	120 mm
Length	450 mm
Depth	450 mm
Refrigerant medium	Ice Dry ice Liquid nitrogen
TBST MAX	93 °C
TBST MIN	-196 °C
Bowls and tubs features	insulating

