



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine light ice tray without lid, vol. 9 l, PVC

Special Price

€113.00 was

€142.00

**plus VAT &
Shipping**

Product Images



Description

The neoLab ice bins are characterized by their light weight and handiness.

They are made of PVC, which makes the containers much lighter than many other ice containers on the market. With a temperature resistance of -196°C to 93°C, the neoLab ice containers are versatile.

The ice containers do not deform and are ideal for use with liquid nitrogen, dry ice, alcohol or salt solutions.

Additional Information

No.	2-4996
Manufacturer (Brand)	neoLabLine
EAN	4058072158316
Transport temperature	Room temperature
Color	Black
Material	Polyvinyl chloride (PVC. non-specific)
Volume	9 l
Wide	365 mm
Width/inside	315 mm
Form	rectangular
Height	120 mm
Height inside	95.0 mm
Length	440 mm
Refrigerant medium	Ice Dry ice Liquid nitrogen
Length inside	380 mm
TBST MAX	93 °C
TBST MIN	-196 °C
Bowls and tubs features	insulating

