



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf IsoPack and IsoRack Set, for 0.5 mL tubes, 0 °C

€58.90
plus VAT &
Shipping

Product Images



Description

IsoPack and IsoRack set, contains one IsoPack and one IsoRack each, for 0.5 mL tubes, 0 °C - IsoTherm system: IsoRack (work rack), IsoSafe (insulated box) and two IsoPacks (cold packs) - Ideal for cooling, transporting and storing frozen samples - Dry incubation principle reduces the risk of contamination often associated with classic water bath incubations for safe sample handling - Holds up to 24 reaction tubes (e. g. E.g. 0.5 mL or 1.5 mL/2.0 mL) for high flexibility - Two versions of IsoPack (cold pack) are available: White holds 0 °C for up to 6 hours. Blue holds -21 °C for up to 3 hours - Racks are stackable, autoclavable and centrifugable in the MTP rotor

Additional Information

No.	VB-0973
Manufacturer (Brand)	Eppendorf
VGKL number	296938860
EAN	4058072411152
Gross weight	0.371
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.33400 kg
Country of origin	Germany
Packing width	0.14000 m
Packing height	0.08000 m
Packing depth	0.11000 m
Packaging volume	2000.000000 ccm
Customs tariff number	39269097

