



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf IsoPack and IsoRack Set, for 1.5/2.0 mL tubes, -21 °C

**€58.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

IsoPack and IsoRack set, contains one IsoPack and one IsoRack each, for 1.5/2.0 mL tubes, -21 °C - IsoTherm system: IsoRack (working rack), IsoSafe (insulating box) and two IsoPacks (cold packs) - Ideal for cooling, transporting and storing frozen samples - Dry incubation principle reduces the risk of contamination often associated with classic water bath incubations for safe sample handling - Holds up to 24 reaction tubes (e. g. E.g. 0.5 mL or 1.5 mL/2.0 mL) for high flexibility - Two versions of IsoPack (cold pack) are available: White holds 0 °C for up to 6 hours. Blue holds -21 °C for up to 3 hours - Racks are stackable, autoclavable and centrifugable in the MTP rotor

## Additional Information

---

No.	VB-0981
Manufacturer (Brand)	Eppendorf
VGKL number	296938874
EAN	4058072411220
Gross weight	0.46
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.41400 kg
Country of origin	Germany
Packing width	0.14000 m
Packing height	0.08000 m
Packing depth	0.11000 m
Packaging volume	2000.000000 ccm
Customs tariff number	39269097

