



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf IsoSafe and IsoPack, contains 1 IsoSafe and 3 IsoPack, for 0.5 mL tubes, 0 °C

**Price on
request
plus VAT &
Shipping**

Product Images



Description

IsoSafe and IsoPack, contains 1 IsoSafe and 3 IsoPack, for 0.5 mL tubes, 0 °C - IsoTherm system: IsoRack (working rack), IsoSafe (insulating box) and two IsoPacks (cold packs) - Ideal for cooling, transporting and storing frozen samples - Dry incubation principle reduces the contamination risk often associated with classic water bath incubations for safe sample handling - Holds up to 24 reaction tubes (e.g. E.g. 0.5 mL or 1.5 mL/2.0 mL) for high flexibility - Two versions of IsoPack (cold pack) are available: White holds 0 °C for up to 6 hours. Blue holds -21 °C for up to 3 hours - Racks are stackable, autoclavable and centrifugable in the MTP rotor

Additional Information

No.	VB-2743
Manufacturer (Brand)	Eppendorf
VGKL number	296938820
EAN	4058072672317
Gross weight	1.2
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.08000 kg
Country of origin	Germany
Customs tariff number	39269097

