



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf IsoSafe and IsoPack, for 1.5/2.0 mL tubes, 0 °C

**€195.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

IsoSafe and IsoPack, contains 1 IsoSafe and 3 IsoPack, for 1.5/2.0 mL tubes, 0 °C - IsoTherm system: IsoRack (work rack), IsoSafe (insulating box) and two IsoPacks (cold packs) - Ideal for cooling, transporting and storing frozen samples - Dry incubation principle reduces the risk of contamination often associated with classic water bath incubations for safe sample handling - Holds up to 24 reaction tubes (e. g. E.g. 0.5 mL or 1.5 mL/2.0 mL) for high flexibility - Two versions of IsoPack (cold pack) are available: White holds 0 °C for up to 6 hours. Blue holds -21 °C for up to 3 hours - Racks are stackable, autoclavable and centrifugable in the MTP rotor

## Additional Information

---

No.	VB-0977
Manufacturer (Brand)	Eppendorf
VGKL number	296938826
EAN	4058072411190
Gross weight	1.2
Transport temperature	Room temperature
HSNumber	39269097
Storage temperature	Room temperature
Net weight	1.08000 kg
Country of origin	Germany
Packing width	0.22000 m
Packing height	0.19000 m
Packing depth	0.17000 m
Packaging volume	7000.000000 ccm
Customs tariff number	39269097

