



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf IsoTherm-System® Starter Set, for 0.5 mL tubes

€179.00
plus VAT &
Shipping

Product Images



Description

IsoTherm-System® starter set, contains IsoSafe, IsoRack, IsoPack 0 °C and IsoPack -21 °C, for 0.5 mL tubes - IsoTherm-System: IsoRack (working rack), IsoSafe (insulating box) and two IsoPacks (cold packs) - Ideal for cooling, transporting and storing frozen samples - Dry incubation principle reduces the risk of contamination often associated with classic water bath incubations for safe sample handling - Holds up to 24 reaction tubes (e. g. E.g. 0.5 mL or 1.5 mL/2.0 mL) for high flexibility - Two versions of IsoPack (cold pack) are available: White holds 0 °C for up to 6 hours. Blue holds -21 °C for up to 3 hours - Racks are stackable, autoclavable and centrifugable in the MTP rotor

Additional Information

No.	VB-0968
Manufacturer (Brand)	Eppendorf
VGKL number	296938811
EAN	4058072411138
Gross weight	1.2
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.08000 kg
Country of origin	Germany
Packing width	0.22000 m
Packing height	0.19000 m
Packing depth	0.17000 m
Packaging volume	8000.000000 ccm
Customs tariff number	39269097

