



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf IsoTherm System® Starter Set, for 1.5/2.0 mL tubes

€179.00
plus VAT &
Shipping

Product Images



Description

IsoTherm-System® starter set, contains IsoSafe, IsoRack, IsoPack 0 °C and IsoPack -21 °C, for 1.5/2.0 mL tubes - IsoTherm system: IsoRack (working rack), IsoSafe (insulating box) and two IsoPacks (cold packs) - Ideal for cooling, transporting and storing frozen samples - Dry incubation principle reduces the risk of contamination often associated with classic water bath incubations for safe sample handling - Holds up to 24 reaction tubes (e. g. E.g. 0.5 mL or 1.5 mL/2.0 mL) for high flexibility - Two versions of IsoPack (cold pack) are available: White holds 0 °C for up to 6 hours. Blue holds -21 °C for up to 3 hours - Racks are stackable, autoclavable and centrifugable in the MTP rotor

Additional Information

| | |
|-----------------------|------------------|
| No. | VB-0976 |
| Manufacturer (Brand) | Eppendorf |
| VGKL number | 296938818 |
| EAN | 4058072411183 |
| Gross weight | 1.2 |
| HSNumber | 39269097 |
| Storage temperature | Room temperature |
| Transport temperature | Room temperature |
| Net weight | 1.08000 kg |
| Country of origin | Germany |
| Packing width | 0.22000 m |
| Packing height | 0.19000 m |
| Packing depth | 0.17000 m |
| Packaging volume | 8000.000000 ccm |
| Customs tariff number | 39269097 |

