



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf PCR Cooler 0.2 mL, pink

Special Price

**€52.24** was

**€64.50**

**plus VAT &  
Shipping**

### Product Images



## Description

---

PCR Cooler 0.2 mL, pink - Handling system for the preparation, protection, transport and storage of sensitive samples - More safety for your samples - Clearly understandable temperature display: PCR Cooler changes color when temperature exceeds 7 °C - For PCR tubes such as reaction tubes, tube strips or plates for flexibility in tube selection - Dry incubation principle reduces contamination risk for samples - Cools an entire 96-well PCR plate for more than one hour at 0 °C (with two hours of pre-cooling at -20 °C) for protection of your samples

## Additional Information

---

No.	VB-0983
Manufacturer (Brand)	Eppendorf
VGKL number	296910023
EAN	4058072411237
Gross weight	0.3
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.27000 kg
Country of origin	Germany
Packing width	0.14000 m
Packing height	0.08000 m
Packing depth	0.12000 m
Packaging volume	2000.000000 ccm
Customs tariff number	39269097

