



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle laboratory tray, PP white, inside LxW 180x240 mm, 3 l

**€10.95** As low as  
**€9.55**  
plus VAT &  
Shipping

### Product Images



## Description

---

The laboratory trays are practical, stackable all-purpose trays, mechanically resilient, break-proof, torsion-resistant. The trays have a strong, rounded edge. Suitable as a collection tray for environmentally hazardous, aggressive substances, as an experiment tray, etching tray, storage tray for instruments, medical tray, photo developer tray, goods display tray and as a food display tray. PP, white Manufactured by injection molding, not thermoformed Space-saving stackable Tray not suitable for heating cabinet or dishwasher Temperature resistant up to +60 °C (short-term up to +80 °C)

## Additional Information

---

No.	KL-0476
Dual use indicator	false
Manufacturer (Brand)	Bürkle
VGKL number	230571180
EAN	4058072264116
Transport temperature	Room temperature
autoclavable at 121°C	No
Net weight	0.22550 kg
Material	Polypropylene (PP)
Volume	3 l
Country of origin	Italy
sterile	No
Customs tariff number	39269097
Length	250 mm
Depth	310 mm
Depth inside	240 mm
Length inside	180 mm

