




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 8-channel reservoir PP natural, 10 pcs./pack

**€15.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Small-volume, colored reservoir made of polypropylene for 8-channel pipettes.

The reservoirs can be connected to each other.

The color differences make it easier to identify liquids, reagents, samples and rinsing agents.

Dimensions:

90 x 20 x 20 mm (L x W x H)

## Additional Information

---

No.	7-4751
Manufacturer (Brand)	neoLab
EAN	4058072041205
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
sterile	No
suitable for	Single-channel pipettes Multi-channel pipettes
Number of chambers	without chambers
Number of channels pipette	8
Filling quantity	100 ml
Type reservoir	Reservoir tub

