



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand reagent reservoir 60 ml f. multichannel pipettes, single sterile, 100 pcs./pack

€181.80
plus VAT &
Shipping

Product Images



Description

Reagent reservoir made of PP, well translucent, volume 60 ml.

For multichannel pipettes with 8 or 12 channels.

Autoclavable at 121°C.

Without lid.

Individually sterile packed.

Additional Information

No.	BN-1637
Manufacturer (Brand)	Brand
old neoLab article no.	7-4621
EAN	4058072209537
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polypropylene (PP)
sterile	Yes
suitable for	Single-channel pipettes Multi-channel pipettes
Bowls and tubs features	individually packed
Number of chambers	without chambers
Number of channels pipette	12
Type of sterility	radiation sterile
Filling quantity	60 ml
Type reservoir	Reservoir tub

