



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Moonlab® Reagent reservoir 25ml, sterile, white,  
PS, individually packed, 100 pcs**

**€125.10  
plus VAT &  
Shipping**

## Product Images



## Description

---

This white sterile reagent reservoir is well suited for liquid handling in aseptic environments, such as sterile safety cabinets. Important e.g. for cell culture. Available in 25 ml

- Material: polystyrene
- individually packed
- sterilized by gamma radiation
- different volumes

## Additional Information

---

|                            |  |
|----------------------------|--|
| No.                        | 4-0015   |
| Manufacturer (Brand)       | Moonlab  |
| EAN                        | 4058072084936                                  |
| Transport temperature      | Room temperature                               |
| autoclavable at 121°C      | No   |
| Color                      | White  |
| Material                   | Polystyrene (PS)                               |
| sterile                    | Yes  |
| suitable for               | Single-channel pipettes Multi-channel pipettes |
| Bowls and tubs features    | individually packed                            |
| Figure                     | 1  |
| Number of chambers         | without chambers                               |
| Number of channels pipette | 8  |
| Type of sterility          | radiation sterile                              |
| Filling quantity           | 25 ml  |
| Type reservoir             | Reservoir tub                                  |

