



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Moonlab® Reagent reservoir 100ml, sterile,  
white, PS, individually packed, 50 pcs**

**€80.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

This white sterile reagent reservoir is well suited for liquid handling in aseptic environments, such as sterile safety cabinets. Important e.g. for cell culture. Available in 100 ml

- Material: polystyrene
- Packing: individually packed, 50 pcs./pack
- sterilized by gamma radiation
- different volumes

## Additional Information

---

No.	4-0017
Manufacturer (Brand)	Moonlab
EAN	4058072084912
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
sterile	Yes
suitable for	Single-channel pipettes Multi-channel pipettes
Bowls and tubs features	individually packed
Number of chambers	without chambers
Figure	3
Number of channels pipette	8
Type of sterility	radiation sterile
Filling quantity	100 ml
Type reservoir	Reservoir tub

