



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® Reagent Reservoir 55ml, 12 compartments, mixed colors, PVC, 5 x 5 pcs/pck, (clear, blue, green, red, yellow)

**€18.90
plus VAT &
Shipping**

Product Images



Description

Reagent reservoir with 12 channels. The reservoir can be filled with reagent solution on both sides. One side has a basin, the other side has 12 channels. Ideal for multi-channel pipettes.

- Material: PVC
- Volume: 55 ml
- Packaging: in the colors clear, blue, green, red, yellow

Additional Information

No.	4-0018
Manufacturer (Brand)	Moonlab
EAN	4058072084929
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky) blue green red yellow
Material	Polypropylene (PP)
sterile	No
suitable for	Single-channel pipettes Multi-channel pipettes
Bowls and tubs features	can be used on both sides
Number of chambers	12 without chambers
Number of channels pipette	12
Filling quantity	50 ml
Type reservoir	Reservoir tub

