



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® pipetting tray for 8-channel pipettes,
PS non-sterile, 100 pcs./pack**

**€110.00
plus VAT &
Shipping**

Product Images



Description

Disposable reagent reservoirs for multichannel pipettes.

Made of white PS, non-sterile or sterile.

With collection channel at the bottom for reduced reagent loss.

Graduated volume markings on the inner wall.

Stable stand due to reinforced side walls and widened bottom.

Packed non-sterile in bags of 100, sterile in bags of 5.

Additional Information

No.	7-3220
Manufacturer (Brand)	neoLab
EAN	4058072038700
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
sterile	No
suitable for	Single-channel pipettes Multi-channel pipettes
Number of chambers	without chambers
Number of channels pipette	8
Filling quantity	6 ml
Type reservoir	Reservoir bowl

