



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® transport tray 20 l, 565 x 365 x 110 mm

€37.40
plus VAT &
Shipping

Product Images



Description

Transport and storage containers made of gray PP, closed at the bottom and sides.
Stackable as desired . Corresponds to Euronorm dimensions.
Resistant to many acids and alkalis.

Basic dimensions:

outside: 600 x 400 mm,

inside: 565 x 365 mm.

Accessories: Grid insert for transport tray (Order No. 1-7543 ff.)

Lid for transport tray (Order No. 2-1349)

In combination with the grid inserts, the ideal means of transport for steep-breasted sample bottles.

Additional Information

No.	2-1144
Manufacturer (Brand)	neoLab
EAN	4058072055653
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	gray
Material	Polypropylene (PP)
Volume	20 l
Wide	400 mm
Width/inside	365 mm
Form	angular
Height	120 mm
Height inside	110.0 mm
Length	600 mm
Depth	600 mm
Length inside	565 mm
Bowls and tubs features	stackable/assemblable with special chemical resistance with handle grid insert available as accessory
Type tub/ box	Transport tray/box

