




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® insulating gloves natural latex, size 11

Price on
request

plus VAT &
Shipping

Product Images



Description

Finger gloves with high resistance to heat and cold and excellent mechanical strength.

- permanently loadable: -40°C to +90°C
- short-time loadable: -70°C to +120°C

Additional Information

No.	2-4224
Manufacturer (Brand)	neoLab
EAN	4058072002565
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	beige
Material	Latex
sterile	No
Length	400 mm
Features gloves	waterproof with special mechanical resistance with special chemical resistance with special protection against thermal risks with special protection against cold extra strong extra long angeraut CE-marked
TBST MAX	120 °C
TBST MIN	-70 °C
Thick fingers	1 mm
Thick palm	1
Thick cuff	1 mm
Gender (clothes, shoes, gloves)	Men
Size	XXL
Glove serves	Personal security
Glove is resistant to	B: Acetone K: Sodium hydroxide 40% L: Sulfuric acid 96%
Gloves provide protection against	Heat Cold Chemicals (limited)
Glove size Germany (Unisex)	11
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Protective clothing for the protection of	Hands/ arms
Strength/ Weight	moderately difficult
Type/ Shape Glove	5-Finger

resistant against

Gasoline/ Oil|Water

with special protection against cold

Yes

standardized according to DIN standard (DIN EN)

DIN EN 374-1:2003-12|DIN EN 388:2003-12

