



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® chain glove, medium, red ribbon

Special Price

**€65.00** was

**€169.00**

plus VAT &  
Shipping

### Product Images



## Description

---

Short five-finger glove made of stainless steel ring mesh (AISI 316 L) with textile strap for optimum protection against cuts and punctures.

Lightweight and ergonomic, are almost not noticed by the user even in the most difficult working conditions.

The glove is reversible and can therefore be worn on both sides. Easy to clean and disinfect.

The textile strap at joint height has a push button on both sides and can be individually adjusted in width.

## Additional Information

---

No.	2-4298
Manufacturer (Brand)	neoLab
EAN	4058072083472
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver red
Material	Stainless steel
sterile	No
Features gloves	with special mechanical strength with special chemical resistance suitable for foodstuffs latex-free extra strong wearable with both hands breathable antibacterial
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Gender (clothes, shoes, gloves)	Unisex
Size (small, medium, large)	medium
Glove serves	Personal security
Gloves provide protection against	Lacerations stab wounds
Protective clothing for the protection of	Hands/ arms
Strength/ Weight	easy
Type/ Shape Glove	5-Finger
with special protection against cold	No
standardized according to DIN standard (DIN EN)	DIN EN 1082-3:2000-06

