



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® nitrile protective gloves cat. III, velourized, size 9

€9.37
plus VAT &
Shipping

Product Images



Description

Medium protection category III (high risk) gloves.

The gloves provide full chemical protection according to EN 374 and are resistant to solvents, fuels, oils and greases.

High protection against cuts and punctures.

With antibacterial treatment to protect against rashes, eczema and fungal diseases. Also antistatic.

Additional Information

No.	2-4291
Manufacturer (Brand)	neoLab
EAN	4058072003012
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	green
Material	Nitrile
sterile	No
Length	32 mm
Features gloves	waterproof with special mechanical resistance with special chemical resistance with special protection against microorganisms/bacterial contamination suitable for foodstuffs latex-free extra long antistatic antibacterial angeraut CE-geke
Size	L
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Thick fingers	0.4 mm
Thick palm	0.4
Thick cuff	0.4 mm
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Glove is resistant to	J: n-Heptane K: Sodium hydroxide 40% L: Sulfuric acid 96%
Gloves provide protection against	Chemicals (protection for at least 30 minutes) Microorganisms/ bacteriological contamination Cuts Sting wounds
Glove size Germany (Unisex)	9
Protective clothing for the protection of	Hands/ arms
Protection level chemicals	Level 2 (> 30 min)
Strength/ Weight	moderately difficult
Type/ Shape Glove	5-Finger

resistant against	Water gasoline/oil animal fats
with special protection against cold	No
standardized according to DIN standard (DIN EN)	DIN EN 374-1:2003-12 DIN EN 388:2003-12

