




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® cut protection gloves with knobs, size 10

€28.50
plus VAT &
Shipping

Product Images



Description

Medium-weight cut protection gloves made of para-aramid synthetic fiber, an extremely resistant material, with two integrated steel cores.

Heat impact is decisively delayed; short-term heat resistance up to 100°C.

Weight: approx. 140 g/pair

Additional Information

| | |
|---|---|
| No. | 2-4234 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072002633 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | blue yellow |
| Material | Cotton Polyphenylene terephthalamide (PPTA. Kevlar®) Steel Polyvinyl chloride (PVC. unspecified) |
| sterile | No |
| Length | 28 mm |
| Features gloves | with special mechanical strength with nubs latex-free extra strong breathable angeroughened CE-marked |
| Size | XL |
| Protection category (PPE) | II Protective equipment for defense against hazards (not CAT. I and III) |
| Gender (clothes, shoes, gloves) | Men |
| Glove serves | Personal security |
| Gloves provide protection against | Heat cuts stab wounds |
| Glove size Germany (Unisex) | 10 |
| Protective clothing for the protection of | Hands/ arms |
| Strength/ Weight | moderately difficult |
| Type/ Shape Glove | 5-Finger |
| with special protection against cold | No |
| standardized according to DIN standard (DIN EN) | DIN EN 388:2003-12 |

