




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® acid protection gloves black, size 7 1/2 - 8, pair

€10.10
plus VAT &
Shipping

Product Images



Description

Protective gloves made of polychloroprene, resistant to many chemical media (acids, alkalis, aliphath. solvents, organic oils and fats).

- anatomic shape
- antistatic
- velour inside
- slightly ribbed on the gripping surface

Additional Information

No.	2-4301
Manufacturer (Brand)	neoLab
EAN	4058072003067
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Polychloroprene (CR)
sterile	No
Length	30 mm
Features gloves	waterproof with special chemical resistance latex-free extra long antistatic angeraut CE-marked
TBST MAX	120 °C
TBST MIN	-30 °C
Size	M
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Glove is resistant to	J: n-Heptane K: Sodium hydroxide 40% L: Sulfuric acid 96%
Gloves provide protection against	Chemicals (protection for at least 30 minutes)
Glove size Germany (Unisex)	8
Protective clothing for the protection of	Hands/ arms
Protection level chemicals	Level 2 (> 30 min)
Type/ Shape Glove	5-Finger
resistant against	Gasoline/ oil water animal fats
with special protection against cold	No
standardized according to DIN standard (DIN EN)	DIN EN 374-1:2003-12 DIN EN 388:2003-12

