




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® acid protection gloves black, size 8 1/2 - 9, pair

**€10.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Protective gloves made of polychloroprene, resistant to many chemical media (acids, alkalis, aliphath. solvents, organic oils and fats).

- anatomic shape
- antistatic
- velour inside
- slightly ribbed on the gripping surface

## Additional Information

No.	2-4302
Manufacturer (Brand)	neoLab
EAN	4058072003074
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Polychloroprene (CR)
sterile	No
Length	30 mm
Features gloves	waterproof with special chemical resistance latex-free extra long antistatic angeraut CE-marked
Gender (clothes, shoes, gloves)	Unisex
Size	L
Glove serves	Personal security
Glove is resistant to	J: n-Heptane K: Sodium hydroxide 40% L: Sulfuric acid 96%
Gloves provide protection against	Chemicals (protection for at least 30 minutes)
Glove size Germany (Unisex)	9
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Protective clothing for the protection of	Hands/ arms
Protection level chemicals	Level 2 ( > 30 min)
TBST MAX	120 °C
TBST MIN	-30 °C
Type/ Shape Glove	5-Finger
resistant against	Gasoline/ oil water animal fats
with special protection against cold	No
standardized according to DIN standard (DIN EN)	DIN EN 374-1:2003-12 DIN EN 388:2003-12

