



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® arm warmers made of aramid, up to 200°C, with thumb hole, pair

**€42.00
plus VAT &
Shipping**

Product Images



Description

Arm warmers made of aramid with thumb hole.
Protection against contact heat up to 200°C.
Cut resistant, washable, double knitted.
CE KAT II, EN 407, EN 388.
Length: 35 cm

Additional Information

No.	8-4003
Manufacturer (Brand)	neoLab
EAN	4058072050818
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Aramid (Meta-Aramid)
sterile	No
Length	35 mm
Features gloves	with special protection against thermal risks with thumb hole latex-free extra strong extra long can be worn on both hands angeraut CE-marked
TBST MAX	200 °C
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Gloves provide protection against	Heat
Protective clothing for the protection of	Hands/ arms
resistant against	Water
with special protection against cold	No
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 407:2004-11

