



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® arm cuffs PP/PE laminate, white, length 38 cm

€131.70
plus VAT &
Shipping

Product Images



Description

Arm warmers, 38cm long with elastic on both sides.
Suitable for the food industry, can also be used under laminar flows.
Offer protection against aggressive substances and fine dusts.
Lint-free, abrasion-resistant and breathable.

- Material PP/PE laminate
- water repellent
- tear-resistant

Additional Information

No.	2-2485
Manufacturer (Brand)	neoLab
EAN	4058072154646
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP) Polyethylene (PE)
Length	380 mm
Features gloves	with special chemical resistance
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Gender (clothes, shoes, gloves)	Unisex
Protective clothing for the protection of	Hands/ arms

