



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® hot/cold gloves, carbon fiber, size 9,
long, length approx. 20 cm, pair**

**€100.60
plus VAT &
Shipping**

Product Images



Description

The well-known and well-proven knitted gloves made of carbon-silicon fiber.

- heat resistant up to 95°C, cold resistant up to -80°C
- Five-finger gloves
- forearm length, 36 cm

Ideally suited for work in autoclaves, freezers or with dry ice.

Not suitable for liquid nitrogen.

Complies with EN 407 4/1/4/1.

Long arm cuff

Additional Information

No.	2-4076
Manufacturer (Brand)	neoLab
EAN	4058072002053
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	dark gray
Material	Carbon Silicon Fiber
sterile	No
Features gloves	with special protection against thermal risks with special protection against cold latex-free extra long breathable angeroughened
TBST MAX	95 °C
TBST MIN	-80 °C
Size	L
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Gloves provide protection against	Heat Cold
Glove size Germany (Unisex)	09. May
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	5-Finger
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 407:2004-11

