



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® 5-finger heat protection gloves, Kevlar, up to 350°C, 35 cm long, pair

€58.50
plus VAT &
Shipping

Product Images



Description

Double insulated knitted glove made of Kevlar®, coarse knit.
5 fingers with lining, size 10.
CE Cat. II, EN 407, EN 388
Contact heat up to 350°C
Length: 35 cm

Additional Information

No.	8-4004
Manufacturer (Brand)	neoLab
EAN	4058072050825
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polyphenylene terephthalamide (PPTA, Kevlar®)
sterile	No
Length	35 mm
Features gloves	with special mechanical strength with special protection against thermal risks with special protection against cold latex-free extra long breathable CE-marked
TBST MAX	350 °C
Gender (clothes, shoes, gloves)	Men
Size	XL
Glove serves	Personal security
Gloves provide protection against	Heat
Glove size Germany (Unisex)	10
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	5-Finger
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 407:2004-11

