



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® 5-finger heat protection gloves with
cotton insulation, aramid, up to 350°C, 36 cm
long, pair**

**€88.60
plus VAT &
Shipping**

Product Images



Description

Knitted **or woven gloves made of aramid.**

Cut resistant. Contact heat up to 350°C.

CE KAT II, EN 407, EN 388.

Order No. 8-4008

5-finger glove, knitted, 36 cm long

Additional Information

No.	8-4008
Manufacturer (Brand)	neoLab
EAN	4058072050863
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Aramid (Meta-Aramid)
sterile	No
Length	36 mm
Features gloves	with special mechanical strength with special protection against thermal risks with special protection against cold latex-free extra strong extra long CE-marked
TBST MAX	350 °C
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Gloves provide protection against	Heat
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	5-Finger
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 407:2004-11

