



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® 5-finger heat protection gloves, Preox-Aramid, up to 650°C, 36 cm long, pair

€91.70
plus VAT &
Shipping

Product Images

Description

Heat protective gloves made of Preox-Aramid with cotton lining insulation.

Contact heat up to 350°C.

Cut resistant, good handling.

CE KAT II, EN 407, EN 388.

Additional Information

| | |
|---|---|
| No. | 8-4011 |
| Manufacturer (Brand) | neolab |
| EAN | 4058072084332 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | dark gray |
| Material | Aramid (Para-Aramid). Preox fiber cotton |
| sterile | No |
| Length | 36 mm |
| Features gloves | with special mechanical strength with special protection against thermal risks with special protection against cold latex-free extra strong extra long angeraut CE-marked |
| TBST MAX | 650 °C |
| Protection category (PPE) | II Protective equipment for defense against hazards (not CAT. I and III) |
| Gender (clothes, shoes, gloves) | Unisex |
| Glove serves | Personal security |
| Gloves provide protection against | Heat |
| Protective clothing for the protection of | Hands/ arms |
| Type/ Shape Glove | 5-Finger |
| with special protection against cold | Yes |
| standardized according to DIN standard (DIN EN) | DIN EN 388:2003-12 DIN EN 407:2004-11 |

