



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® high temperature glove, five-finger, up to 650°C, pair

**€120.80
plus VAT &
Shipping**

Product Images



Description

Gloves for the highest demands and temperatures!

Gloves for use especially where it is particularly hot, e.g. for heating and muffle furnaces (protection cat. III).

Five-finger gloves, 300 mm long

- made of preoxaramide fabric
- Five-finger gloves, cut-resistant, with insulation
- short for temperatures up to 650°C
- Length: approx. 300 mm

Additional Information

No.	2-4390
Manufacturer (Brand)	neoLab
EAN	4058072003272
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	dark gray
Material	Aramid (para-aramid). Preox fiber
sterile	No
Length	30 mm
Features gloves	with special mechanical strength with special protection against thermal risks with special protection against cold latex-free extra strong extra long long-skinned
TBST MAX	650 °C
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Gloves provide protection against	Heat
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	5-Finger
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 407:2004-11

