



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® high temperature protective gloves,
mitten, up to 900°C, pair**

**€174.50
plus VAT &
Shipping**

Product Images



Description

Gloves for highest demands and temperatures!

Gloves for use especially where it is particularly hot, e.g. for heating and muffle furnaces (protection cat. III).

Mittens with thumb, 300 mm long

- made of high temperature fabric
- multi-layer insulation, fabric does not stick to glass
- short for temperatures up to 900°C
- Length: approx. 300 mm

Additional Information

No.	2-4391
Manufacturer (Brand)	neoLab
EAN	4058072003289
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Aramid (meta-aramid) fiberglass
sterile	No
Length	30 mm
Features gloves	with special mechanical strength with special protection against thermal risks with special protection against cold latex-free extra strong extra long can be worn ambidextrously angry
TBST MAX	900 °C
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Gloves provide protection against	Heat
Protection category (PPE)	III Complex protective equipment for high risks. to protect against lethal hazards or serious irreversible damage to health
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	Mitten
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 407:2004-11

