



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Hot Hand Heat Protection

€30.80
plus VAT &
Shipping

Product Images



Description

The thick-walled silicone rubber is heat-resistant up to approx. 250°C, but also retains its flexibility down to -57°C.

The two grip openings for thumb and fingers fit any hand size.

The grip surface has knobs and is therefore non-slip.

With eyelet for hanging.

Additional Information

No.	2-4098
Manufacturer (Brand)	neoLab
EAN	4058072002084
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Silicone
sterile	No
Features gloves	waterproof with special mechanical strength with special protection against thermal risks with special protection against cold with nubs suitable for food latex-free extra strong can be worn ambidextrously
TBST MAX	250 °C
TBST MIN	-57 °C
Protection category (PPE)	I Simple protective equipment for low hazards
Gender (clothes, shoes, gloves)	Unisex
Glove serves	Personal security
Gloves provide protection against	Heat Cold
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	Mitten
with special protection against cold	Yes

