



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryo gloves, waterproof, wrist sound, size XL

€280.00
plus VAT &
Shipping

Product Images



Description

100% waterproof cryogenic gloves for work in liquid nitrogen gas phase, argon, helium, oxygen, liquid gas on closed systems, cryogenic refrigerators etc.

Also for work in hot atmosphere or autoclave.

- Protection from -160° to +150°C
- suitable for handling infectious or toxic substances
- with seamless inner glove,
PTFE-coated
- autoclavable

Not suitable for immersion in liquid nitrogen!

Conform to EN 511, 420 and 388.

Additional Information

No.	2-4261
Manufacturer (Brand)	neoLab
EAN	4058072002893
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Cotton Polytetrafluoroethylene (PTFE)
sterile	No
Length	30 mm
Features gloves	waterproof with special mechanical resistance with special chemical resistance with special protection against thermal risks with special protection against cold latex-free suitable for work with liquid nitrogen extra strong extra long
TBST MAX	150 °C
TBST MIN	-160 °C
Size	XL
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Gender (clothes, shoes, gloves)	Men
Glove serves	Personal security
Gloves provide protection against	Heat Cold Cuts Sting Wounds
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	5-Finger
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 420:2010-03 DIN EN 511:2006-07

