




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Cryo gloves, waterproof, 40 cm long, size XL

**€265.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

100% waterproof cryogenic gloves for work in liquid nitrogen gas phase, argon, helium, oxygen, liquid gas on closed systems, cryogenic refrigerators etc.

Also for work in hot atmosphere or autoclave.

- Protection from -160° to +150°C
- suitable for handling infectious or toxic substances
- with seamless inner glove,  
PTFE-coated
- Autoclavable

Not suitable for immersion in liquid nitrogen!

Conform to EN 511, 420 and 388.

## Additional Information

No.	2-4213
Manufacturer (Brand)	neoLab
EAN	4058072002497
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Cotton Polytetrafluoroethylene (PTFE)
sterile	No
Length	40 mm
Features gloves	waterproof with special mechanical resistance with special chemical resistance with special protection against thermal risks with special protection against cold latex-free suitable for work with liquid nitrogen extra strong extra long
TBST MAX	150 °C
TBST MIN	-160 °C
Gender (clothes, shoes, gloves)	Men
Size	XL
Glove serves	Personal security
Gloves provide protection against	Heat Cold Cuts Sting Wounds
Protection category (PPE)	II Protective equipment for defense against hazards (not CAT. I and III)
Protective clothing for the protection of	Hands/ arms
Type/ Shape Glove	5-Finger
with special protection against cold	Yes
standardized according to DIN standard (DIN EN)	DIN EN 388:2003-12 DIN EN 511:2006-07 DIN EN 420:2010-03

